

## Analytical Industries Inc.

2855 Metropolitan Place, Pomona, CA 91767 USA 🔸 Tel: 909-392-6900, Fax: 909-392-6900, e-mail: sales-medical@aii1.com,Website: www.aii1.com

Oxygen Sensor Model	PSR-11-915-2i	
Measuring Range	0-100% O2	
Signal Output (1)	45-65 mV	
Response Time 90%	9 Seconds	[28.96]
Accuracy Full Scale (2)	<u>+</u> 1%	Ø1.140
Accuracy FS Over Operating Range (3)	<u>+</u> 5%	
Drift % Signal Month (2)	< 1%	
Linearity (2)	<u>+</u> 1%	
Humidity Non-Condensing	0-99% RH	[27.18]
Temperature Coefficient	Compensated	Ø1.070
Operating Temperature	0-45° C	[21,97]
Operating Pressure	Constant 5-30 PSIG	[17.53]
Flow Rate	.1-10 LPM	NOTE: FOR OPTIMUM PERFORMANCE TORS SHIP EACH C DOWN OF
Expected Life (1)	60 Months	[1.91] [28.14] NOTE: FOR OPTIMUM PERFORMANCE HIS END FACING DOWN OR HORIZONTAL
Storage (4)	0 to 40°C	
Warranty (5)	18 Months	
Electrical Connections	3 Ring PCB Gold Plated	
Packaging	Polyetylene bag in 3" x 3" x 1.5" paper box	
Recommended Storage	< 6 months	

- (1) In air @ 25° C and 1 atm. Expected life varies inversely with changes in O2 concentration, temperature, pressure.
- (2) Under constant conditions. Calibrate as near sample conditions as possible. Calibrate with 95-100% oxygen for measurements > 40% oxygen.
- (3) Once the signal output stabilizes following a step change in temperature.
- (4) Recommended 0-40° C. Sensors may be stored at -10 to 55° C on an intermittent basis only.
- (5) Ex-factory unless agreed to in writing. Under normal application conditions sensor is warranted to be free of defects in material and workmanship for the period specified above proved the sensor is properly installed and operated. Analytical Industries Inc. shall not be liable for buyer's negligence, misapplication, alteration or abuse. The sole remedy for a sensor determined to be defective by Analytical Industries Inc. is limited to replacing the defective sensor.